



THUNDER Phospho SMAD3 (S423/S425) + Total SMAD3 TR-FRET Cell Signaling Assay Kit

CATALOG NUMBERS KIT-SMAD3PT-500 (500 tests)
400 points for phospho-SMAD3
and 100 points for total SMAD3

Store at -80°C
For research use only.
Not for use in diagnostic procedures.

[Visit Product Page](#)

PRODUCT DESCRIPTION

This assay kit measures intracellular levels of **phospho-SMAD3 (S423/S425)** and **total SMAD3** protein in cell lysates using a simple, rapid and sensitive immunoassay based on the homogeneous (no-wash) THUNDER™ TR-FRET technology. The kit is compatible with both adherent and suspension cells.

SPECIFICITY

This assay kit contains two specific and selective antibody pairs, one that recognizes **SMAD3** phosphorylated at **Ser423** and **Ser425**, and another that recognizes **total** (both phosphorylated and unphosphorylated) **SMAD3**.

SPECIES REACTIVITY

Human and Mouse (Swiss-Prot Acc.: P84022; Entrez-Gene Id: 4088).

Other species should be tested on a case-by-case basis.

TR-FRET ASSAY PRINCIPLE

The **Phospho-SMAD3 (S423/S425) + Total SMAD3** assay kit is a homogeneous time-resolved Förster resonance energy transfer (TR-FRET) sandwich immunoassay (Figure 1). The THUNDER™ Cell Signaling assay workflow consists of 3 steps (Figure 2). Following cell treatment, cells are first lysed with the specific Lysis Buffer provided in the kit. Then **Phospho-SMAD3 (S423/S425) + Total SMAD3** in the cell lysates are detected in separate wells with two pairs of fluorophore-labeled antibodies in a simple "add-incubate-measure" format (single-step reagent addition; no wash steps). For detection of the phosphorylated protein, one antibody is labeled with a donor fluorophore (Europium chelate; Eu-Ab1) and the second with a far-red acceptor fluorophore (FR-Ab2). The same approach is used for the second antibody pair detecting the total protein (Eu-Ab3 and FR-Ab4). The binding of the two matched labeled antibodies to distinct epitopes on the target protein (either **phospho-SMAD3** or **total SMAD3**) takes place in solution and brings the two dyes into close proximity. Excitation of the donor Europium chelate molecules with a flash lamp (320 or 340 nm) or a laser (337 nm) triggers a FRET from the donor to the acceptor molecules, which in turn emit a TR-FRET signal at 665 nm. Residual energy from the Eu chelate generates light at 615 nm. The signal at 665 nm is proportional to the concentration of **Phospho-SMAD3 (S423/S425)** and **Total SMAD3** in the cell lysate. Data can be expressed as either the signal at 665 nm or the 665 nm/615 nm ratio.




Figure 1 Schematic representation of the TR-FRET cell signaling assay principle.




Figure 2 Assay workflow using the 2-plate (transfer) protocol.

KIT COMPONENTS

	500 points*
Eu-labeled phospho-SMAD3 (S423/S425) antibody (Eu-Ab1)	20 µL
Acceptor-labeled phospho-SMAD3 (S423/S425) antibody (FR-Ab2)	80 µL
Eu-labeled total-SMAD3 antibody (Eu-Ab3)	5 µL
Acceptor-labeled total-SMAD3 antibody (FR-Ab4)	20 µL
Lysis Buffer 1 (5X)	5 mL
Detection Buffer (10X)	250 µL
Positive control cell lysate	200 µL
Phosphatase Inhibitor Cocktail (100X)	250 µL

* The number of assay points is based on an assay volume of 20 µL in half-area 96-well or low-volume 384-well assay plates using the kit components at the recommended concentrations (refer to the User Manual).

TECHNICAL DATA SHEET

Phospho-SMAD3 (S423/S425) + Total SMAD3


VALIDATION DATA

This assay kit has been validated for the relative quantification of phospho SMAD3 (S423/S425) and total SMAD3 in A549 cell lysates using the 2-plate assay protocol.

- Adherent cells were cultured overnight in a 96-well tissue culture plate (RPMI +10% FBS).
- Following cell treatment, the media was removed and cells were lysed with the 1X **Lysis Buffer 1** (50 μ L) supplemented with the 100X Phosphatase Inhibitor Cocktail diluted at 1X.
- Following a **30-min** incubation at room temperature (RT) on an orbital shaker (400 rpm), lysates (15 μ L) were transferred to a 384-


- well assay plate followed by addition of the labeled antibodies Eu-Ab1 and FR-Ab2 (5 μ L) for detection of phospho SMAD3 (S423/S425) or Eu-Ab3 and FR-Ab4 (5 μ L) for detection of total SMAD3.
- The plate was incubated at RT for **18 hours** and the TR-FRET signal was recorded at 665 and 615 nm (PHERASTAR® FSX; laser excitation).

STIMULATION OF PHOSPHO-SMAD3 (S423/S425) IN A549 CELLS




A549 cells (75,000 cells/well; in triplicate) were incubated with serial dilutions of TGF-β1 for 15 min at 37°C. Data show that treatment of A549 cells with TGF-β1 stimulates phosphorylation of SMAD3 at S423/S425, but does not have an effect on the levels of total SMAD3.

INHIBITION OF PHOSPHO-SMAD3 (S423/S425) IN A549 CELLS




A549 cells (75,000 cells/well; in triplicate) were incubated with serial dilutions of the inhibitor SB 431542 for 15 min at 37°C. Cells were then stimulated with 0.5 nM of TGF-β1 for 15 min at 37°C. Data show that treatment of A549 cells with SB 431542 inhibits phosphorylation of SMAD3 at S423/S425 by TGF-β1 but does not have an effect on the levels of total SMAD3.

A549 CONTROL LYSATE TITRATION (QC TEST) PHOSPHO-SMAD3 (S423/S425)



A549 CONTROL LYSATE TITRATION (QC TEST) TOTAL SMAD3



Quality Control: the total SMAD3 assay kit is routinely tested against TGF-β1 treated A549 lysates. A549 cells were cultured in a T175 flask to 85% confluence and stimulated with 3 nM of TGF-β1 for 15 min at 37°C. Following cell lysis using 4 mL of 1X Lysis Buffer 1, lysates were serially diluted with 1X Lysis Buffer 1 and tested in triplicate. Data show a linear relationship between lysate dilutions and TR-FRET ratio values.



bioauxilium

BETTER TOOLS. REAL DISCOVERIES.

FOR MORE INFORMATION ON DEVELOPING AND OPTIMIZING TR-FRET CELL SIGNALING ASSAYS, CONSULT THE USER MANUAL.

7171 Frederick-Banting, Suite 2115, Montréal (Québec) H4S 1Z9 CANADA bioauxilium.com 1-866-610-6933 x 103